

i-Optik

NEW



TCS-760

Ophthalmic chair and stand

Features

- Well designed and layout, Smooth movements and good practical functions.
- Operated by control panel
- Electric up/down controls for both chair and table
- Chair arm can be moved up and down, and can be rotated 90° in the highest
- Sliding board on the desktop can hold optical instrument
- Phoropter arm can attach manual vision tester or computerized vision tester.
- Trial lens set can be placed inside of the table drawer

Specifications

Table	Table size	1060×494×30mm
	Elevation range for table	698—848mm
	Sliding board size	290×420mm
	Input voltage	AC220V/50Hz, 110V/60Hz
	Weight of table	150Kg
Chair	Elevation range for chair	520-670mm
	Weight of chair	45Kg

